

BATHING CHAIRS

BATHING CHAIRS



Foldable bathing chair

Frame made from white coloured aluminium with ergonomic seat with cut-away for better personal hygiene.

Item nr. 302010



Foldable bathing chair with WC

Frame made from white coloured aluminium with a removable toilet bucket with lid and carrying handle.

Item nr. 302011

COMMUNE CHAIRS

COMMUNE CHAIRS

The commode and shower chair has a frame made of white anodized aluminium and comes as standard with 4 pieces of 5" castors with parking brakes. The seat, backrest and armrests are made of soft black PU. The foot-rests are adjustable in height and are removable. Comes complete with toilet bucket with lid and carrying handle.

The commode and shower chair can be equipped with 24" rear wheels with parking brakes. In the 24" wheel kit are included 2 brakes with clamps and brackets for fitting of the wheels on the frame. Pre-drilled holes makes the mounting of the wheels very easy.

Anti-tip is optional.



Technical data	Small 302015	Standard 302018	Large 302019	XL 302020	Tilt 302024	Tilt L 302025	Diff. by 24"
Seat width	38 cm	45 cm	50 cm	55 cm	45 cm	50 cm	---
Seat depth	38 cm	43 cm	43 cm	43 cm	43 cm	43 cm	---
Seat height (adjust.)	50 - 55 cm	50 - 55 cm	50 - 55 cm	50 - 55 cm	55 - 60 cm	55 - 60 cm	---
Pass over height	42 - 47 cm	42 - 47 cm	42 - 47 cm	42 - 47 cm	---	---	---
Pass over width	32 cm	40 cm	45 cm	50 cm	---	---	---
Total width	49 cm	56 cm	61 cm	66 cm	56 cm	61 cm	+ 8 cm *
Total depth	85 cm	93 cm	93 cm	93 cm	110-120 cm	110-120 cm	+ 5 cm
Total height (adjust.)	96-101 cm	96-101 cm	96-101 cm	96-101 cm	97-141 cm	97-141 cm	---
Max. load	130 Kg	150 Kg	150 Kg	200 kg	130 Kg	130 Kg	---
Weight	11,00 Kg	12,90 Kg	13,50 Kg	14,10 kg	15,10 kg	15,70 kg	+ 2,90 kg
Tilt angle	---	---	---	---	-5° to+25°	-5° to+25°	---

* Total width: Difference for Kakadu Tilt (302024 or 302025) with 24" wheels mounted: + 14 cm



BATHING PRODUCTS



Dealer

SHOWER CHAIR

The Mobilex shower chairs all have a frame made of powder-coated high-quality aluminium with a white fixed plastic seat. A soft seat cover is available as accessory for the triangular chairs (picture B1 and B2) The round shower stool in picture C1 has a rotating seat pad, the model in picture C2 has a white fixed plastic seat.



Description	Pic.	Dimensions	Max. load
Wall mounted shower chair	A	44 x 54,5 x 43,5-51 cm	100 kg
Triangular shower chair	B1	50 x 45 cm	120 kg
Triangular shower chair	B2	44 - 46 cm	150 kg
Round shower chair (Ø x H)	C1	33 cm x 35,5 - 53,5 cm	130 kg
Round shower chair (Ø x H)	C2	33 cm x 40-50 cm	130 kg

SHOWER CHAIR

The shower chair in picture D has a grey seat made of soft PUR and is available with white or silver-grey frame and in various sizes. The legs of all models are height adjustable.



Description	Pic.	Dimensions	Max.load
Shower chair	D	52,5-62,5 x 53 x 60-77 cm	150 kg
Budget shower chair	E	30 x 50 x 37,5-50,5 cm	115 kg

BATH TUBE SEATS

The Mobilex bath tub seat (F) and the swivel bather (G) both have a hygienic cut-out. The frame made of silver-coloured aluminium. Seat and backrest are made of grey soft PUR. Back and armrests are easily removable (clips). The backrest for the bath tube seat (F) is optional.



Description	Pic.	Dimensions	Max.load
Bath tube seat	F	76-91 x 42 x 23 cm	150 kg
Swivel bather	G	71 x 54 x 57 cm	130 kg

SLIDING BATH TUBE SEAT



The sliding bathtub seat has an aluminium frame with white plastic seat, back and soap holder. The seat is height-adjustable and mounted on roller bearings and can be easily slid over the bathtub wall with the help of the handle.

Sliding bathtub seat	Dimensions
Width x depth	107-111 x 49-53 cm
Seat width x seat depth	50 x 39,5 cm
Drive-over bathtub edge	45,5 - 55,5 cm
Weight	6,5 kg
Max. load	130 kg

BATH BOARDS



The bath boards are made of injection moulded and recyclable white plastic. They enable the user to sit comfortably on the bathtub. A smooth surface makes cleaning easier.

Description	Pic.	Dimensions	Max.load
Bath board	H1	71 x 32 x 19 cm	190 kg
Bath board	H2	70 x 35 x 17 cm	150 kg
Bath board XL	--	76 x 35 x 18 cm	150 kg