

The "**Kudu**" lightweight rollator has a black or silver coloured aluminium frame that can be folded sideways with in folded state reduces the required space to a minimum. The rollator is equipped with a back strap and comes also with a handy removable bag made of tear-resistant nylon. The bag has carrying handles and can also be used as a shopping bag. Easy to expand and collapse the rollator with its large 200 mm wheels is suitable for use both outdoors as indoors. The comfortable seat is made of tear-resistant nylon. The push handles can be adjusted in height and allow a comfortable and secure progress. The rollator comes with two lockable brakes on the rear wheels, with tip-assist, a crutch holder and a bike bell.

Mobilex rollator Kudu	Standard Item no. 312075	Medium Item no. 312076	Small Item no. 312077
<i>Height of the push handles (adjustable)</i>	80 - 98,5 cm	79 - 88 cm	69 - 80 cm
<i>Seat size (H x W x D)</i>	65 x 42 x 20,5 cm	57 x 42 x 20,5 cm	50 x 42 x 20,5 cm
<i>Size of the folded rollator (H x W x D)</i>	82 x 29 x 67 cm	79 x 27 x 68 cm	67 x 27 x 72 cm
<i>Total width</i>	61 cm	61 cm	61 cm
<i>Inner width between push handles</i>	45,5 cm	45,5 cm	45,5 cm
<i>Length (depth) opened</i>	67 cm	68 cm	70 cm
<i>Weight</i>	8,0 kg	7,9 kg	7,7 kg
<i>Max. load (rollator / bag)</i>	150 kg / 5 kg	150 kg / 5 kg	150 kg / 5 kg
<i>Brakes</i>	2 lockable with parking function		
<i>Wheel sizes</i>	Front and back 200 x 35 mm		



RK1a. Mobilex "Kudu" Rollator



RK1b. The folded "Kudu" Rollator

Standard equipment



RK1c. Bike bell



RK1d. Tray



RK1e. Cane-holder



RK1f. Tip-assist