

Indoor rollator



317410



317411



317412



317414

New, updated model:

- Lighter weight
- New wheel design
- Aluminium frame which folds easily for storage
- Lightweight rollator that is very stable and easy to manoeuvre
- Equipped with four smooth-running wheels
- Equipped with two lockable brakes, which are operated by means of a handle
- The rollator is tested and approved for a user weight of up to 100 kg
- The rollator is only approved for indoor use



317416

As standard the Jaguar LIGHT comes with:

- a bag
- a removable tray

Jaguar LIGHT	Item no. 317410-12+14+16-17 Various colours	Item no. 317415 SUPREME
Height of push handles	83,5 - 94,5 cm	90,3 - 96 cm
Folded rollator size (H x W x D)	83,5 x 56 x 26 cm	90,3 x 56 x 26 cm
Width	56 cm	56 cm
Length opened / folded	66 cm / 26 cm	66 cm / 26 cm
Weight	5,84 kg	0
Max. load (rollator / bag)	100 kg / 5 kg	100 kg / 5 kg
Brakes	2 lockable with parking function	
Wheel sizes	Front and back: 7" (175 x 27 mm)	
Turning circle	80,8 cm	
Tube diameter	Handle tube Ø21*1,5*402mm, Front connecting tube: Ø19*1*398mm, Fixed tube rear frame 19*1*395mm, Frame front Ø 18,5x36,3/frame back Ø26.5x39,3mm	
Material	Aluminium frame, plastic (PA+GF), basket of fabric	



317417



317415 - SUPREME