

## *Verticalizer and walking frame in modern design*



*Accessory - item no. 313380:  
An adjustable supporting seat belt  
can be ordered as an accessory.*

**Bison is a walking aid for people with walking difficulties or balance problems. The walking frame is also highly suitable for rehabilitation.**

The walker has various functions and adjustment options:

- The two soft polyurethane forearm supports can be easily adjusted in width to the user and help the user maintain an upright position while walking, by resting the entire forearm on the forearm supports.
- The height of the walker is easily adjusted up or down using a gas spring.
- The handles with brake arms can be adjusted in both angle and depth.
- The rear wheels have a directional brake and the front wheels have a total brake.
- All adjustments are easily made without the use of tools.
- The narrow width of the walker makes it possible to get through narrow passages.
- The walker is made of powder-coated aluminium and is intended for indoor use only.

Bison	313300
<i>Width of the frame at the rear wheels</i>	66,0 cm
<i>Width of arm rests min. - max.</i>	66,0 - 80,0 cm
<i>Distance between handles</i>	34,0 cm
<i>Height min. - max.</i>	102,2 - 143,0 cm
<i>Length</i>	86,0 cm
<i>Weight</i>	18,3 kg
<i>Max. load</i>	150 kg

Scan the QR code  
and learn more:

