

ELEPHANT AIR



The Elephant Air rollator with large pneumatic tyres is one of the best outdoor rollators currently available. This rollator is of high quality and has extra large tyres for more stability on difficult surfaces.

The outdoor rollator weighs only 10.68 kg, has a back strap, a walking stick holder and a closed bag. With its large pneumatic tyres, the Elephant Air rollator is suitable for both indoor and outdoor use and is tested and approved for a user weight of up to 150 kg.



The rollator can also be folded very easily with a flick of the wrist by pulling the seat upwards on a loop. Then you can also fold the wheels so that the rollator is even smaller and can be easily transported.

Rollator Elephant Air	Crème Royal 317001	Matte black 317002
Height of the push handles (adjustable)	87,3 - 104,8 cm	
Recommended user heights	150 - 195 cm	
Seat height	62 cm	
Total width opened / folded	67 cm / 28,0 cm	
Inner width between push handles	45,5 - 40,7 cm	
Length (depth)	80,5 cm	
Weight (+ bag 380 g)	10,68 kg	
Max. load (rollator / bag)	150 kg / 5 kg	
Brakes	2 lockable with parking function	
Wheel sizes	Front and back: 12" / 291 mm	