

NEW ARRIVAL



Absolutely unique position on the market

PANDA COMFORT

Item no. 312601

- Dark blue aluminium frame that folds easily for storage and transport.
- Very stable and easy to manoeuvre thanks to the smooth-running wheels.
- The large soft wheels provide good shock absorption.
- In driving position it is secured against unwanted unlocking
- With smooth-running and lockable brakes.



- The optimal construction of the frame and the handles is the basis for the best possible operating comfort.
- A wide range of accessories available.
- Tested and approved for a user weight of up to 150 kg.

As standard the Panda Comfort comes with:

- a transport lock
- a removable, closed shopping bag
- a back belt
- a crutch holder

Scan the QR code
and learn more:

