

# NEW ARRIVAL



## SEAL

The Seal is a simple and solid wheelchair that makes the use comfortable and reliable without increasing the complexity. The wheelchair has the following features:

- Foldable frame
- Raisable/ lowerable
- Armrests
- Removable footrests
- Parking brakes
- Quick release axles (option for 24" wheel version)



### Configuration options for the wheelchair:

- Width: 40, 45 and 50 cm
- Wheel: 24" or 12" PU tires
- Axle: Quick release axles for the 24" wheel version available

Scan the QR code  
to learn more:



***Only available in containers (approx. 130 pcs.) - delivery time 4 months***  
***From 79€ FOB depending on the constellation***

# SEAL

Available with **customized logos** and in **selected colors**.  
Containers available from **MOQ of 130 pieces**:



With your  
logo here



## Spare parts and accessories for the wheelchair:

- Brakes
- 24" rear wheel made of PU with handrim
- 8" x 1" front wheel made of PVC
- Front fork
- Padded seat cover
- Legrest
- Adjustable leg rest (optional)
- Padded back cover
- Footrest
- Footplate made of black plastic
- Armrest with side panel
- Armrest PU
- Safety belt (optional)
- Tire pressure brake (optional)
- Anti-tip (optional)