

# NEW ARRIVAL



## *Absolutely unique position on the market*

# PANDA ROYAL

Item no. 312602

## Carbon fibre outdoor rollator

- Purple/black carbon fibre frame
- Total weight of only 8,4 kg
- The rollator folds easily for storage and transport
- Very stable and easy to manoeuvre thanks to the smooth-running wheels
- The large soft wheels provide good shock absorption
- In driving position it is secured against unwanted unlocking
- With smooth-running and lockable brakes



- The optimal construction of the frame and the handles is the basis for the best possible operating comfort
- A wide range of accessories available
- Tested and approved for a user weight up to 150 kg

**As standard the Panda Royal comes with:**

- a transport lock
- a removable, closed shopping bag
- a back belt
- a crutch holder

Scan the QR code  
and learn more:



***Expected on stock Mid-November 2024***

# NEW ARRIVAL



*Absolutely unique position on the market*

## PANDA ROYAL

Item no. 312602

### Carbon fibre outdoor rollator



Panda Royal	312602
<i>Height of the push handles (adjustable)</i>	80,8 - 99,3 cm
<i>Recommended user heights</i>	160 - 200 cm
<i>Seat height</i>	59 cm
<i>Seat dimensions (width x depth)</i>	45,5 x 19 cm
<i>Inner width between push handles</i>	45,2 - 46,2 cm
<i>Length (depth)</i>	71,7 cm
<i>Weight</i>	8,4 kg
<i>Max. load (rollator / bag)</i>	150 kg / 5 kg
<i>Brakes</i>	2 lockable with parking function
<i>Wheel sizes</i>	Front and back: 283 mm x 42/43 mm