

TURTLE

Item no. 271045 / 271145

Turtle is a simple, yet solid and lightweight transport wheelchair made in a stylish silver aluminium frame. Turtle has just enough features to make the usage comfortable and reliable without making it too complicated.

Turtle has the following features:

- Foldable frame
- Armrests
- Removable footrests
- Removable, washable seat cover
- Parking and attendant brakes
- Anti-tip pedal
- Seat width 45 cm
- Tested and approved for a user weight up to 150 kg



Item no. 271045



Item no. 271145



Available in two different models:

- With 6" front castor wheels and 14" soft PU rear wheels attached
- With 6" front castor wheels and 22" PU rear wheels with quick-release axles attached

Scan the QR code
to learn more:



— **Expected on stock in the beginning of December 2024** —

NEW ARRIVAL



TURTLE

Item no. 271045 / 271145



271045



271145

Turtle	271045	271145
Seat-width	45 cm	45 cm
Seat-depth adjustable	42-48 cm	
Seat height adjust. min./max.	46,5 cm – 49,0 cm	
Total width	62,4 cm	66,3 cm
Backrest inclination	10°	
Weight	12,5 kg	16 kg
Max. user weight	150 kg	