

# NEW ARRIVAL



## Comfort and flexibility in one chair

# TOILET CHAIR

Item no. 302013

- Can be used both as a:
  - Toilet chair
  - Raised toilet seat
  - Shower chair
- Made of lightweight material
- With soft, removable PU seat and back
- With folding armrests for easy access in and out of the chair
- Can easily be adjusted in height without the use of any tools



As standard the chair comes with:

- Toilet bucket with lid
- Small side tray\*
- Toilet paper holder\*

\*can be mounted on both sides



Technical Data	Toilet chair with toilet bucket 302013
Total size (LxWxH)	51,6-58,7 x 64,5-67,5 x 72-92 cm
Folded size (LxWxH)	51,5 x 64,5 x 72 cm
Width between armrests	44 cm
Height adjustable	72-92 cm
Seat height adjustable	44,5 - 64,5 cm
Weight	7,0 kg
Max. load	150 kg

Scan the QR code  
and learn more:



**Expected on stock in early July 2025**