

NEW ARRIVAL



XXL model of the popular and solid:

KAKADU

Item no. 302021

Commode and shower chair



- Frame made of white anodized aluminum
- Comes with 4 pieces of 5" twin castors with parking brakes
- The seat, backrest and armrests are made of soft black PU
- The footrests are easily adjustable in height and are removable
- Comes with a toilet bucket with carrying handle, a lid and a black PU seat cover
- Max. load: 270 kg

*Please note
24" rear wheels cannot be fitted to the Kakadu XXL,
as the chair is not approved for this.*



Scan the QR code
and learn more

NEW ARRIVAL



XXL model of the popular and solid:

KAKADU

Item no. 302021

Commode and shower chair



Technical data	Small 302015	Standard 302018	Large 302019	XL 302020	Difference with 24"	XXL 302021
Seat width	38 cm	45 cm	50 cm	55 cm	---	61 cm
Seat depth	38 cm	43 cm	43 cm	43 cm	---	51 cm
Seat height adjust.	50 - 55 cm	50 - 55 cm	50 - 55 cm	50 - 55 cm	---	51-53 cm
Seat opening	15,5 x 21 cm	19 x 26 cm	19 x 26 cm	19 x 26 cm	---	19 x 26 cm
Pass over height **	42,5 - 47,5 cm	43 - 48 cm	42 - 47 cm	40 - 45 cm	---	41-43 cm
Pass over width	32 cm	40 cm	45 cm	50 cm	---	51 cm
Seat/Back angle	90 ° / 12°				---	90 ° / 12°
Distance between armrests	39,5 cm	46,5 cm	51,5 cm	56,6 cm	---	62,5 cm
Backrest height	46 cm				---	46 cm
Total width	49 cm	56 cm	61 cm	66 cm	+ 9,5 cm	72,5 cm
Total depth	85 cm	93 cm	93 cm	93 cm	+ 19 cm	93 cm
Total height adjust.	96-101 cm	96-101 cm	96-101 cm	96-101 cm	---	99-101 cm
Max. load	130 Kg	150 Kg	150 Kg	200 kg	---	270 kg
Weight	11,00 Kg	12,90 Kg	13,50 Kg	14,10 kg	+ 4,20 kg	19,75 kg
** +5 cm with the new accessory 30212M						